



Non Adjustable Splint

- Custom-fitted to individual's measurements for maximum stability
- Ideal for consistent joint size or strong reinforcement needs
- Distribute pressure evenly and prevent shifting
- Accurate initial sizing is necessary



Adjustable Splint

- Allow for widening or tightening due to fluctuations in joint size
- Suitable for patients with EDS, arthritis, or variable swelling
- Quick, at-home adjustments without replacement
- Over time, may provide slightly less stability, and might require re-shaping of the band

ADJUSTABLE SPLINTS = best for flexibility, swelling management, and trial fittings

NON-ADJUSTABLE SPLINTS = best for long-term stability, durability, and precise support